In vitro diagnostic. For self-testing.

IMPORTANT: Please read this insert and your Owner’s Booklet before using. The Owner’s Booklet has full instructions for blood glucose and control solution testing. **Do Not** use new test strips if the vial is open or damaged in any way. This could lead to error messages or inaccurate blood glucose results. Contact LifeScan Customer Service right away at 1 888 567-3003 (English) (7 days a week, 8 a.m. – 10 p.m. Eastern Time. For assistance outside of these hours, please contact your health care professional) if you have questions about your test strips, test results or these instructions.

IMPORTANT SAFETY INSTRUCTIONS

- This meter and lancing device are intended to be used by a single person and should not be shared, even with a family member.
 - After use and exposure to blood, all parts of this kit are considered biohazardous. A used kit can potentially transmit infectious diseases even after you have performed cleaning and disinfection.
- Refer to your Owner’s Booklet for important information about cleaning and disinfecting your meter, lancing device, and caps.

Intended Use

The OneTouch® Verio™ Test Strips are used with the OneTouch® Verio™ Family of Blood Glucose Monitoring Systems (OneTouch® Verio™IQ) for the quantitative measurement of glucose (sugar) in fresh capillary whole blood samples drawn from the fingertips. The OneTouch® Verio™ Family of Blood Glucose Monitoring Systems is intended to be used by a single patient and should not be used for testing multiple patients.

Refer to the Owner’s Booklet that came with your meter for the full system Intended Use.

Storage and Handling

- Store the test strip vial in a cool dry place between 41°F and 86°F. Keep away from direct sunlight and heat. Exposure to temperatures and/or humidity outside the storage limits may cause inaccurate readings.
- The expiration date is printed on the vial.
- The discard date is 6 months after you first open a test strip vial. When you open a new vial, write the discard date on the label.
- **Do Not** use after the expiration or discard date, whichever comes first.
- Store your test strips in the **vial only**. To avoid damage or contamination, **Do Not** put test strips in any other container.
- Open the vial **only** when taking out a test strip for use.
- Close the vial lid tightly right away after removing a test strip. Use each test strip **as soon as you take it out** of the vial.
- **Do Not** use test strips from a vial that is damaged or left open to air.
- **Do Not** get dirt, food, or liquids on the test strip. You may touch the test strip anywhere on its surface with clean, dry hands.
- **Do Not** bend, cut, or alter test strips in any way.
- Test strips are for single use only.
- **Never reuse a test strip that already has blood or control solution on it.**
- Make sure your meter and test strips are about the same temperature before you test.
- Apply only OneTouch® Verio™ Control Solution or a blood sample to the test strip.
- **Do Not** put the used test strip back in the vial after taking a test.
- Used test strips may be considered biohazardous waste in your area. Be sure to follow your health care professional’s advice or local laws when throwing them away.

⚠ WARNING: Keep the test strip vial away from children. Test strips are a choking hazard. Do Not swallow test strips. The test strip vial may contain drying agents that are harmful if inhaled or swallowed and may cause skin or eye irritation. Do Not ingest or swallow any items.

Testing Your Blood Glucose

Read your system Owner’s Booklet for instructions on running a blood test and getting a blood sample.

Test Results

The OneTouch® Verio™ Family of Meters shows results between 20 and 600 mg/dL.

Read these cautions when your test results are higher or lower than what you expect.

⚠ CAUTION:

Low glucose results
If your test result is below 70 mg/dL or a low glucose warning, it may mean low blood glucose (hypoglycemia). You may need to treat low glucose right away. Follow your health care professional’s advice. The result could be due to a test error, but it is safer to treat first, then test again.

Dehydration and low glucose results
You may get false low glucose results if you are severely dehydrated. If you think this is the case, contact your health care professional right away.

High glucose results
If your test result is above 180 mg/dL, it may mean high blood glucose (hyperglycemia). You should think about re-testing. Talk to your health care professional if you are concerned about high blood glucose.

A high glucose warning is shown when your blood glucose level is over 600 mg/dL. You may have very high blood glucose (severe hyperglycemia). Re-test your blood glucose. If the result is a high glucose warning again, this indicates a severe problem with your blood glucose control. Obtain and follow instructions from your health care professional right away.

Repeated unexpected glucose results
If you get repeated unexpected results:

- Run a control solution test to check your system.
- Make sure you follow all instructions in your Owner’s Booklet.
- If your symptoms don’t match your blood glucose results, call your health care professional.

Never ignore symptoms or make significant changes to your diabetes management program without speaking to your health care professional.

Range of Expected Values

- Blood glucose management requires the help of a health care professional. Talk to your health care professional about:
- The target blood glucose range that is right for you.
 - When to test.
 - What your blood glucose results mean.

Expected blood glucose levels for non-pregnant people without diabetes¹:

Time	Range, mg/dL
Fasting	Less than 100
2 hours after meals	Less than 140

Checking the System

Use only OneTouch® Verio™ Control Solutions

- A control solution test is used to check that:
- The meter and test strips are working together properly.
 - You are doing the test correctly.

- Do a control solution test**
- When you open a new vial of test strips.
 - If you think the meter or test strips are not working properly.
 - If you get repeated unexpected blood glucose results.
 - If you drop or damage the meter.

Read your system Owner’s Booklet for instructions on how to run a control solution test.

Limitations of Procedure

- OneTouch® Verio™ Test Strips give accurate results under the following conditions:
- Use only fresh capillary whole blood. **Do Not** use serum or plasma.
 - Hematocrit is the percentage of red blood cells in the blood. Extremes in hematocrit may affect test results. Hematocrit levels below 20% may cause false high readings. Hematocrit levels over 60% may cause false low readings. Talk to your health care professional if you don’t know your hematocrit level.
 - OneTouch® Verio™ Test Strips may be used at altitudes up to 10,000 feet with no effect on test results.
 - Fresh capillary blood samples may be collected into sodium heparin, lithium heparin, potassium EDTA and sodium citrate containing test tubes if the blood is used within 10 minutes. **Do Not** use sodium fluoride/oxalate or other anticoagulants or preservatives.
 - Interferences: Acetaminophen, ascorbic acid (vitamin C), and uric acid (when occurring in blood at high normal or at high therapeutic concentrations) do not significantly affect results. However, abnormally high concentrations in blood may cause inaccurately low results. Tolazamide levels up to 7.5 mg/dL (3x the level found at high therapeutic levels) do not affect the blood glucose results.
 - Test results may be falsely low if the patient is severely dehydrated, in shock, or in a hyperosmolar state (with or without ketosis).
 - Critically ill patients should not be tested with blood glucose meters.
 - Lipemic samples: Cholesterol levels up to 700 mg/dL and triglycerides up to 3000 mg/dL do not affect the results. Grossly lipemic patient samples have not been evaluated and are not recommended for testing with the OneTouch® Verio™ Family of Meters.

⚠ CAUTION: The OneTouch® Verio™ Family of Blood Glucose Monitoring Systems should not be used for patients within 24 hours of receiving a D-xylose absorption test as it may cause inaccurately high results.

Alternate Site Limitations

Refer to the Owner’s Booklet that came with your meter for the full system Intended Use. Some meter systems are not indicated for Alternate Site Testing. Results from forearm or palm samples may differ from results using fingertip samples. This is due to normal differences in blood flow. It does not reflect system accuracy. Be sure to talk to your health care professional before using your forearm or palm for sampling. Differences in blood glucose readings between the fingertip and other sites may occur after eating, insulin medication or exercise. This is because blood sugar levels change faster in the fingertip.

Alternate site testing (AST) should not be used to calibrate continuous glucose monitors (CGMs) or to calculate an insulin dose. AST testing should only be done during steady-state times (when glucose is not changing rapidly).

Please note that the OneTouch® Verio™IQ meter is not indicated for Alternate Site Testing.

Test Principle

The OneTouch® Verio™ Family of Meters is plasma-calibrated for easy comparison of results with lab methods. Glucose in the blood sample mixes with special chemicals on the test strip to produce an electrical current. The strength of the current changes with the amount of glucose in the blood sample. The OneTouch® Verio™ Family of Meters measures the current, calculates your blood glucose level and shows your result.

Reagent Composition

Each test strip contains: flavin adenine dinucleotide dependent glucose dehydrogenase or FAD-GDH (from *Aspergillus* sp.) - 2 U; potassium ferricyanide 41 µg; and other ingredients (buffer, etc.). The vial contains a drying agent.

Performance Characteristics

The performance of the OneTouch® Verio™ Test Strip has been evaluated both in laboratory and in clinical tests.

Measurement Range: The measurement range of the OneTouch® Verio™ Family of Meters is 20 to 600 mg/dL.

System Accuracy

Diabetes experts have suggested that glucose meters should agree within 15 mg/dL of a laboratory method when the glucose concentration is lower than 75 mg/dL, and within 20% of a laboratory method when the glucose concentration is 75 mg/dL or higher. Samples from 100 patients were tested using both the OneTouch® Verio™IQ System and the YSI 2300 Glucose Analyzer laboratory instrument.

System Accuracy Results for Glucose Concentrations <75 mg/dL		
Percent (and number) of meter results that match the laboratory test		
Within ±5 mg/dL	Within ±10 mg/dL	Within ±15 mg/dL
88.2% (45/51)	100.0% (51/51)	100.0% (51/51)

System Accuracy Results for Glucose Concentrations ≥75 mg/dL			
Percent (and number) of meter results that match the laboratory test			
Within ±5%	Within ±10%	Within ±15%	Within ±20%
71.1% (177/249)	94.8% (236/249)	98.0% (244/249)	100% (249/249)

System Accuracy Results across the entire Glucose Range	
Percent (and number) of meter results that match the laboratory test	
Within ±15 mg/dL or ±20%	
100% (300/300)	

Therefore, 100% of the total results obtained with the OneTouch® Verio™IQ System achieved the goal suggested by the diabetes experts.

Regression Statistics			
Samples were tested in duplicate on each of three test strip lots. Results indicate that the OneTouch® Verio™IQ System compares well with a laboratory method.			
# of Subjects	# of Tests	Slope	Intercept (mg/dL)
100	300	0.95	4.81

95% CI Slope	95% CI Intercept (mg/dL)	Std. Error (S _{y.x}) (mg/dL)	R²
0.94 to 0.96	2.76 to 6.87	9.73	0.99

Precision			
Within Run Precision (300 Venous Blood Tests)			
Target Glucose (mg/dL)	Mean Glucose (mg/dL)	Standard Deviation (mg/dL)	Coefficient of Variation (%)
40	42.39	1.21	2.86
100	97.43	1.84	1.89
130	130.24	2.49	1.91
200	196.15	3.96	2.02
350	330.34	7.39	2.24

Results show that the greatest variability observed between test strips when tested with blood is 2.86% or less

Total Precision (600 Control Solution Tests)			
Glucose Level Ranges (mg/dL)	Mean Glucose (mg/dL)	Standard Deviation (mg/dL)	Coefficient of Variation (%)
Low (38-62)	50.31	1.77	3.52
Normal (102-138)	117.14	2.70	2.31
High (298-403)	342.34	8.68	2.53

Data generated using the OneTouch® Verio™IQ meter. The OneTouch® Verio™IQ meter is representative of the OneTouch® Verio™ Family of Meters. Refer to your system Owner’s Booklet for performance data specific to your system.

IMPORTANT: For instructions, technical data, and a list of symbols used, read your system Owner’s Booklet.

References

1. American Diabetes Association, Position Statement, Diagnosis and Classification of Diabetes Mellitus, Diabetes Care 31:555-560, 2008.

OUR COMMITMENT TO YOU

Our goal is to provide you with quality health care products and dedicated customer service. If you are not fully satisfied with this product, or if you have questions about the use of any LifeScan product, contact us 7 days a week, 8 a.m. – 10 p.m. Eastern Time at 1 888 567-3003 (English). If you cannot reach Customer Service, contact your health care professional for advice.

Covered by one or more of the following U.S. patents: 6,179,979, 6,284,125, 6,716,577, 6,749,887, 6,863,801, 7,045,046, 7,498,132, and 7,846,312. Use of these test strips and associated monitoring device is protected under the following U.S. patents: 6,413,410, 7,749,371. Purchase of the associated monitoring device does not act to grant a use license under these patents. Such a license is granted only when the associated monitoring device is used with OneTouch® Verio™ Test Strip. No test strip supplier other than LifeScan, Inc. is authorized to grant such a license. The accuracy of results generated with LifeScan meters using test strips manufactured by anyone other than LifeScan has not been evaluated by LifeScan.

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Division of Cilag GmbH